



P-BAG.

P-BAG, a more convenient and hygienic substitute for Bedpan and Bottle. Suitable for women and men.

- **P-BAG**, is a urine collecting plastic bag, equipped with a Unique membrane which prevent leakage and bad odour. UNI-BAG can easily replace all other as so called open urine container As bedpan and bottles etc. Both men and women can use UNI-BAG at different kind of positions and situations May be used in a Sitting, lying down or standing position by women.
- **P-BAG**, has a unique handle for easy grapping, carrying and Can be hanged up and reused several times up to a volume Of 1.5 litres without leakage and bad odour. No odour emission Is precluded by unique check valve system.
- **P-BAG** gives a better working condition for the staff, no risk of Contamination from urine to staff, no more smell and bad odour, No more leaks as before from open bottles and
- **P-BAG**, is fully disposable - made from **Polyethylene** not **PVC** therefore more environmentally friendly. Moulded, jug Especially for women.
- **P-BAG**, reduce a lot of working time for the staff and Thereby decrease the costs for the hospital, institution etc.

Product fact.

Art No.	10 0005,	10 0007	10 0008.
Number of pieces per Carton.	30 Pcs,	100 Pcs	50 Pcs
Number of Cartons per Pallet.	120	100	72 Cartons.
Volume.	1 500 ml.		
Material.	Polyethylene HDPE and LDPE.		

Clean Analytical System AB

Address: Gårdfarivägen 3 532 38 Skara	Telephone: +46 511-163 70	Fax: +46 511-180 82	E-post: info@casab.se	Website: www.casab.se
--	-------------------------------------	-------------------------------	---------------------------------	---------------------------------